

WHITE PAPER

Medical device recall and FDA letter to health care providers

FLEXIBLE UROLOGICAL ENDOSCOPES AND REPROCESSING METHODS

Ambu White Paper - Flexible Urological Endoscopes

April 2022

On Friday April 1, 2022 Karl Storz Endoscopy–America Inc. initiated an urgent medical device recall notice on certain Karl Storz flexible endoscopes for urological use. Shortly after that the Food and Drug Administration (FDA) issued a letter to health care providers ensuring that all relevant stakeholders were informed about the details of the urgent recall and understood their options.

Prior to the recall, supplemental validation testing of the efficacy of the manual high-level disinfection (HLD) process was performed showing that the required efficacy level of disinfection was not achieved. As a result, the manufacturer has now removed manual and automated HLD as well as liquid chemical sterilization (for specific models) from the Instruction for Use (IFU) of the affected endoscopes.

The FDA has provided recommendations to the end users of the affected devices^{1,2}:

1. Do not use high-level disinfection methods or liquid chemical sterilization to reprocess affected urological endoscopes.
2. Sterilize affected urological endoscopes after each use by using sterilization methods recommended in the instructions for use specific to each device.
3. Do not use affected urological endoscopes if you do not have access to an appropriate sterilization method recommended in the instructions for use.
4. Discuss the benefits and risks associated with procedures involving reprocessed urological endoscopes with your patients.

As early as April 2021, FDA communicated about infections associated with reprocessed urological endoscopes providing recommendations for how to reprocess reusable urological endoscopes³. Based on over 450 Medical Device Reports (MDRs), which described patient infections post procedure or other possible contamination issues associated with reprocessing, FDA raised awareness about proper reprocessing processes according to the manufacturers' IFUs.

WHAT DOES THIS CHANGE IN REPROCESSING METHODS MEAN FOR THE USER AND FOR PATIENT SAFETY?

CURRENT PRACTICE AND UNMET NEEDS

In the United States, high-level disinfection (HDL) is the minimum recommended practice for the reprocessing of reusable flexible cystoscopes⁴. It includes seven separate steps associated with organizational and financial burden: Pre-cleaning, leak testing, manual cleaning, visual inspection, disinfection, storage, and documentation^{4,5}. The overall cost of reprocessing one flexible endoscope with HDL without additional sterilization has been estimated to range between \$114 and \$280⁶.

Sterilization of endoscopes is intended to add an additional safety margin relative to the current reprocessing methodology⁷. However, sterilization of endoscopes is an expensive and time-consuming process, which requires specific equipment, facilities, and trained staff⁸. Furthermore, some evidence suggests that sterilization of endoscopes may not solve the issues of infection and cross-contamination^{9,10}.

IMPLICATIONS OF STERILIZATION AS ONLY POSSIBLE REPROCESSING METHOD

Besides liquid chemical sterilization, other methods of sterilizing endoscopes include the traditional steam sterilization method via autoclave, hydrogen peroxide gas plasma sterilization, and Ethylene Oxide (EtO) sterilization. These methods all have different benefits and limitations. Since most reusable flexible cystoscopes are heat and moisture-sensitive, the traditional method of physical sterilization using steam is not a relevant solution to implement in current practices⁴.

The EtO sterilization method was highlighted by both the FDA and the Environmental Protection Agency (EPA) to have a significant negative impact on the environment and human health. EPA has classified EtO as a hazardous air pollutant subject to tight restrictions within the Clean Air Act, and in general, EtO sterilization is under increasing scrutiny and pressure from both the FDA and the EPA^{11,12}.

In July 2019, FDA announced two innovation challenges aimed at developing new device sterilization methods and reducing emissions connected to EtO sterilization¹². During the same year, the EPA issued a state EPA order to close a facility sterilizing medical products due to unacceptable levels of EtO in the air around the facility¹³.

In addition to the environmental burden, the EtO sterilization method has a long turnaround time of approximately 2 hours and 30 minutes compared to 1 hour and 15 minutes for the conventional HLD method^{6,14}. EtO sterilization has also been found to add an additional per procedure cost of \$339 (2015 USD)¹⁵.

We expect both health systems and manufacturers to continue to be challenged by the FDA and EPA to move away from EtO sterilization to solutions that are less harmful for the environment and human health.

Furthermore, the use of gases and chemicals to sterilize endoscopes are known to cause significant damage to the endoscopes. This imparts additional costs and leads to endoscopes being unavailable due to repairs^{10,15,16}.

Since all recommended sterilization methods may result in an increase in the turnaround time and higher repair rates, it must

be assumed to negatively affect current reusable flexible cystoscopy capacity. Hence, it is expected that health care systems will experience additional costs associated with ensuring a bigger fleet of cystoscopes to perform the same amount of procedures, establishing new reprocessing procedures and training the staff accordingly.

HOW CAN A SINGLE-USE SOLUTION OFFER SUPPORT?

Sterilization of reusable urological endoscopes is cumbersome, costly and time-consuming, and it contributes to the deterioration of the endoscope. What, then, could be an alternative that would still permit adequate patient treatment?

The implementation of a single-use cystoscopes is a solution that can allow health care providers to provide an economically viable, safe and timely option for patient treatments and can be used in facilities where proper sterilization is not applicable.

The cost of implementing single-use cystoscopy differs between facilities. However, the per procedure cost of a cystoscopy procedure has been found to range from \$155-495, which makes single-use cystoscopes an economically sustainable solution¹⁷⁻¹⁹. The clinical performance of single-use cystoscopes has been found to be superior or comparable to reusable cystoscopes^{20,21}. What is more, the continuous availability of single-use cystoscopes makes it possible to maintain your current procedure volume without adding additional resources even if you are affected by the recall, and it helps reduce the risk of procedure delays^{19,22,23}.

In conclusion, single-use cystoscopy is regarded as the preferred alternative to reusable flexible cystoscopy. It provides a safe and simple option for immediate patient treatment on an ongoing basis and/or when reusable endoscopes are unavailable for whatever reason. Single-use cystoscopes offer a superior or comparable clinical outcome, are cost transparent and enable institutions to retain maximum cystoscopy capacity. Due to potential workflow benefits and the elimination of reprocessing and repairs, single-use cystoscopes can free up resources to focus on patient care and other important tasks.

REFERENCES

1. KARL STORZ. Urgent Medical Device Recall Notice Certain KARL STORZ Flexible Endoscopes for Urological Use. https://www.karlstorz.com/cps/rde/xbcr/karlstorz_assets/ASSETS/3680244.pdf (2022).
2. The U.S. Food and Drug Administration (FDA). UPDATE: Change in Reprocessing Methods with Certain Karl Storz Urological Endoscopes - Letter to Health Care Providers. <https://www.fda.gov/medical-devices/letters-health-care-providers/update-change-reprocessing-methods-certain-karl-storz-urological-endoscopes-letter-health-care> (2022).
3. The U.S. Food and Drug Administration (FDA). Infections Associated with Reprocessed Urological Endoscopes - Letter to Health Care Providers. <https://www.fda.gov/medical-devices/letters-health-care-providers/infections-associated-reprocessed-urological-endoscopes-letter-health-care-providers> (2021).
4. American Urological Association. Joint AUA/SUNA White Paper on Reprocessing of Flexible Cystoscopes. <https://www.auanet.org/guidelines/guidelines/flexible-cystoscopes> (2018).
5. Centers for Disease Control and Prevention (CDC) & The Department of Health and Human Services (HHS). Essential Elements of a Reprocessing Program for Flexible Endoscopes - Recommendations of the Healthcare Infection Control Practices Advisory Committee. <https://www.cdc.gov/hicpac/pdf/flexible-endoscope-reprocessing.pdf> (2017).
6. Ofstead, C. L., Quick, M. R., Eiland, J. E. & Adams, S. J. A Glimpse at the True Cost of Reprocessing Endoscopes: Results of a Pilot Project. *International Association of Healthcare Central Service Material* 62-78 (2017).
7. Rutala, W. A., Kanamori, H., Sickbert-Bennett, E. E. & Weber, D. J. What's new in reprocessing endoscopes: Are we going to ensure "the needs of the patient come first" by shifting from disinfection to sterilization? *Am J Infect Control* 47S, A62-A66 (2019).
8. Scotland, K. et al. Evaluation of a Single-Use Flexible Cystoscope: A Multi-Institutional International Study. *J Endourol* 34, 981-986 (2020).
9. Larsen, S. et al. Rate and impact of duodenoscope contamination: A systematic review and meta-analysis. *EClinicalMedicine* 25, (2020).
10. Ofstead, C. L. et al. The effectiveness of sterilization for flexible ureteroscopes: A real-world study. *Am J Infect Control* (2017) doi:10.1016/j.ajic.2017.03.016.
11. United States Environmental Protection Agency (EPA). What EPA Is Doing to Address Ethylene Oxide and to Learn More About the Chemical | US EPA. <https://www.epa.gov/hazardous-air-pollutants-ethylene-oxide/what-epa-doing-address-ethylene-oxide-and-learn-more-about#regulations>.
12. The U.S. Food and Drug Administration (FDA). Ethylene Oxide Sterilization for Medical Devices. <https://www.fda.gov/medical-devices/general-hospital-devices-and-supplies/ethylene-oxide-sterilization-medical-devices#epa> (2022).
13. The U.S. Food and Drug Administration (FDA). Statement from FDA Commissioner Scott Gottlieb, M.D., on steps the Agency is taking to prevent potential medical device shortages and ensure safe and effective sterilization amid shutdown of a large contract sterilization facility. <https://www.fda.gov/news-events/press-announcements/statement-fda-commissioner-scott-gottlieb-md-steps-agency-taking-prevent-potential-medical-device> (2019).
14. Centers for Disease Control and Prevention (CDC). Ethylene Oxide Sterilization. <https://www.cdc.gov/infectioncontrol/guidelines/disinfection/sterilization/ethylene-oxide.html> (2008).
15. Almario, C. v et al. Cost-Utility of Competing Strategies to Prevent Endoscopic Transmission of Carbapenem-Resistant Enterobacteriaceae. *Am J Gastroenterol* 110, 1666-1674 (2015).
16. Omidbakhsh, N., Manohar, S., Vu, R. & Nowruzi, K. Flexible gastrointestinal endoscope processing challenges, current issues and future perspectives. *J Hosp Infect* 110, 133-138 (2021).
17. Su, Z. T., Huang, M. M., Matlaga, B. R., Hutfless, S. & Koo, K. A micro-costing analysis of outpatient flexible cystoscopy: implications for adoption of single-use flexible cystoscopes. *World J Urol* 39, 4275-4281 (2021).
18. Wong, A., Phan, Y. C., Thursby, H. & Mahmalji, W. The First UK Experience with Single-use Disposable Flexible Cystoscopes: An In-depth Cost Analysis, Service Delivery and Patient Satisfaction Rate with Ambu® aScopeTM 4 Cysto. *J Endolum Endourol* 4, e29-e44 (2021).
19. Pietropaolo, A., Hughes, T., Tear, L. & Somani, B. K. Comparison of ureteric stent removal procedures using reusable and single-use flexible cystoscopes following ureteroscopy and lasertripsy. *Cent European J Urol* 73, 342-348 (2020).
20. Whelan, P., Kim, C., Tabib, C., Preminger, G. M. & Lipkin, M. E. Evolution of Single-Use Urologic Endoscopy: Benchtop and Initial Clinical Assessment of a New Single-Use Flexible Cystoscope. *J Endourol* 36, 13-21 (2022).
21. Rindorf, D. K., Besser, V. & Blædel, T. Clinical Performance of the Single-Use Cystoscope aScopeTM 4 Cysto: An evaluation based on initial perceptions from urologist worldwide. *Ambu A/S* (2021).
22. Doherty, D. T., Moinuddin, Z., Grey, B. R. & van Dellen, D. IsirisTM for Ureteric Stent Removal in Renal Transplantation: An Initial Single-Centre Experience of 150 Cases. *Surg Innov* 28, 366-370 (2021).
23. Baston, E. L. et al. Office-based ureteric stent removal is achievable, improves clinical flexibility and quality of care, whilst also keeping surgeons close to their patients. *Cent European J Urol* 71, 196-201 (2018).