

Ambu® Disposable Face Mask

DATASHEET

Description

The Ambu® Disposable Face Mask is an anatomically shaped face mask designed for patient oxygenation and ventilation.

The face mask comes in 6 different sizes from neonate to large adult, and is a single patient multiple use face mask, to minimizing risk of cross contamination.

Features

The mask is available with and without different features such as check valve and hook ring. Further it has a thumb rest, designed to improve grip during handling.

The anatomically shaped cushion allows a tight fit to the patient's face, enabling a tight seal.



The Ambu® Disposable Face Mask connected to the Ambu® SPUR® II Resuscitator, in various sizes.



KEY BENEFITS

Transparent dome allows for patient observation during use

Anatomically shaped flexible cushion, allows a tight fit to the patient's face

Optional with and without hook ring

MR safe

Phthalate-free material

Ambu

ORDERING SPECIFICATIONS WITHOUT HOOK RING

Item no.	Type
000252051	Ambu Disposable Face Mask Neonate
000252052	Ambu Disposable Face Mask Infant
000252053H	Ambu Disposable Face Mask Toddler
000252054H	Ambu Disposable Face Mask Small adult/Child
000252055H	Ambu Disposable Face Mask Medium adult
000252056H	Ambu Disposable Face Mask Large adult
000252951	Ambu Disposable Face Mask w. check valve Neonate
000252952	Ambu Disposable Face Mask w. check valve Infant
000252953H	Ambu Disposable Face Mask w. check valve Toddler
000252954H	Ambu Disposable Face Mask w. check valve Small adult/Child
000252955H	Ambu Disposable Face Mask w. check valve Medium adult
000252956H	Ambu Disposable Face Mask w. check valve Large adult

WITH HOOK RING

Item no.	Type
000252053	Ambu Disposable Face Mask Toddler
000252054	Ambu Disposable Face Mask Small adult/Child
000252055	Ambu Disposable Face Mask Medium adult
000252056	Ambu Disposable Face Mask Large adult
000252953	Ambu Disposable Face Mask w. check valve Toddler
000252954	Ambu Disposable Face Mask w. check valve Small adult/Child
000252955	Ambu Disposable Face Mask w. check valve Medium adult
000252956	Ambu Disposable Face Mask w. check valve Large adult

SIZE

Type	Weight
Neonate	9 g
Infant	14 g
Toddler	25 g
Small adult/Child	33 g
Medium adult	47 g
Large adult	49 g

MATERIALS

Product	
Dome	Polyvinylchloride
Cushion	Polyvinylchloride
Check valve	Polycarbonate
Packaging	
Pouch	PP (Polypropylene)
Box	Cardboard
Environment	
Phthalate-free material	Yes
MR safe	Yes



Ambu® Disposable Face Mask with check valve.



Ambu® Disposable Face Mask without check valve.

Ambu

Ambu A/S

Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
ambu.com