


# Ambu® Head Wedge®



## Key Benefits

- Disposable
- Waterproof and very durable
- Radiolucent ensuring that the patient is effectively immobilized even during examination
- Stores completely flat, saving valuable storage space, it can even be stored attached to the backboard
- Fully adjustable to fit all sizes from children to adults
- Fits all known backboards
- Strong, self-adhesive head strap with medical grade glue ensures a secure hold without any discomfort to the patient

## Ambu® Head Wedge®


The unique Ambu Head Wedge is designed to immobilize the patient's head instantly and effectively when placed on a backboard.

It is applied easily and quickly in just four simple steps. With its innovative construction it enables a perfect fit to the patient in one go, saving time at the scene when seconds count.

See interactive flash on our home page [www.ambu.com](http://www.ambu.com)

## Specifications

Description	
Accommodates head sizes	50 cm to 63 cm in circumference (19.7" to 24.8")
Dimensions	(L x W x T) 400 x 165 x 12 mm (15.8" x 6.5" x 0.5")
Weight	100 g (3.5 ounces)
Recommended operating temperatures	-5°C - +40°C (-23°F - +104°F)
Storage temperatures	-40°C - +70°C (-40°F - +158°F)
Adjustable to head sizes in	30 increments
Color	Yellow
Environmental properties	Waterproof and UV-resistant
Head comfort	Foam pad
Tape liners	Easy pull tape liner
Headstrap	Medical grade adhesive



Ambu® Head Wedge®

## Material


Product	Material
Head wedge base	Polypropylene
Actuator	Polyamid (PA) (Nylon)
Foam pad	Polyethylene (Recyclable PE)
Headstrap	Polyester (UP)




## Ordering specifications

Item no.	Description
000 264 034	Ambu Head Wedge (with one strap), 10/pack
000 264 036	Ambu Head Wedge (with two straps), 10/pack
000 264 035	Ambu Head Wedge (no strap), 10/pack

Available packaged with Perfit ACE and Mini Perfit ACE



Ambu® Head Wedge®



Ambu® Perfit ACE®  
Ambu® Mini Perfit ACE®